

Disinfectible PC-Keyboard for Purity and Hygiene

AK-C8100 Series Medical Key*



Highlights:

- Hygiene Keyboard for Clinical Wipe Disinfection
- Standard PC Key Layout with 104/105 Keys
- 360° Safeguard against Dirt and Splash Water
- Flat Keys for Safe Routine Disinfection
- Contained Key Membrane Made of Silicone
- Seamless Surface for Complete Wipe Disinfection
- Easy Key Actuation for Perfect Typing
- Effective Usability and Ergonomics
- Minimum Life Cycle Cost by Efficient Data Entry
- Exchange Membrane for Machine Reprocessing
- Protection Class IP68 and RF-Wireless as Option



Features & Benefits:

The AK-C8100 is a high quality, high performance, and highly reliable PC-keyboard. The purification and disinfection capacity makes it a perfect fit in hospitals, at practitioners and in microbiological laboratories. The silicone membrane protects keyboard and key field against ingress of dirt and fluids. The AK-C8100 impresses with congenial key actuation and precise key guidance. With its standard PC key layout, the AK-C8100 is highly recommendable for heavy duty type writing in hygiene critical applications.

The AK-C8100 increases protection against infection and limits risk and liability of the health care provider:

- Thorough wipe disinfection of key field and keyboard reduces microbial burden in the hospital
- Silicone key membrane with low key profile and contained, seamless surface offers effective protection against cross contamination, proliferation of microbes and formation of biofilms
 - ✓ No vector for nosocomial pathogens by building an effective germ barrier at the workplace
 - ✓ No microbial growth on debris and residues from disinfectants in hollow spaces
 - ✓ No spread of germs and surface recontamination by leakage of polluted cleaning suds
- CleanFunction allows to switch off the keys while cleaning ("Fn"+"Clean on/off")
- Approved for clinical wipe disinfection according to VAH-List (Germany)
- Machine reprocessing of exchange membrane at 93°C and 134°C (in addition to wipe disinfection)
or: Full sealing IP68, respectively sealed membrane IP65 (RF-Wireless)

Specification (Wired)	
Compatibility	All Operating Systems that support Standard USB resp. PS/2 Keyboard
Temperature Range	Storage -20° to +70° C, Operation 0° to +50° C
Keyboard Protection Membrane	Exchangeable Silicone Membrane (Standard) or Fully Sealed IP68 (Option)
Chemical Resistance	Wipe Disinfection according to Guidelines Cleaning and Disinfection DGHM/VAH
Interface	USB 1.1 (1.2 MBit/s)
Supply Voltage	5 V +5%/-10% at USB Port resp. 5V +/-5% at PS/2 Port
Current Consumption	Approx. 35 mA (USB: Low Power Device)
Connector	USB Type A resp. PS/2 Mini-DIN
Dimensions	Approx. 449 x 142 x 26 mm
LED Status Indicators	Caps Lock, Num Lock, Scroll Lock, CleanFunction
Key Switch Technology	Scissors Balancing Key Guidance with Contained Contact System
Key Travel	Approx. 2,5 +/- 0,5 mm
Actuation Force	Approx. 60 grams (similar to Notebook style Keyboards)
Agency Approvals, Certificates	CE, RoHS, EN-60601-1-2

*No Medical Device according to 93/42/EWG

Part# / Keyboard:		AK-C8100F-zzz-y/xx (Wired)	AK-C8100F-Fzz-y/xx (RF)
Part# / Exchange Membrane:		AK-C8100-y/xx	Specification see Datasheet AK-C8100F-FU
Identifier zzz = Version/Option U1 = Standard (USB), PS/2 upon request UVS = IP68 Fully Sealed (Option) FU1 = RF, FUS = RF Sealed (Option)	Identifier y = Color W = White (Standard) WB = White Housing/ Black Membrane	Identifier xx = Country Version US = US English BE = Belgium CH = Swiss IT = Italian GE = German FR = French NOR = Nordic UK = UK English SP = Spanish Other Versions upon Request	

Medical Key

Exchangeable Membrane for Routine Disinfection & Sterilization



Work



Sterilize



Exchange

© 2016 Active Key GmbH & Co. KG. Errors and omissions excepted. Warranties excluded.

**Keyboards for Industry
Commerce and Service**



Active Key


Active Key GmbH & Co. KG
Brunnenäcker 6 91257 Pegnitz, Germany

Phone: +49 (0)9241 48337 0
Fax: +49 (0)9241 48337 29

www.activekey.eu
www.activekeyboards.com

info@activekey.de

Clean & Write



AK-C8100.

P/N (US English, White):

AK-C8100F-U1-W/US (Standard Wipe Disinfection)

AK-C8100F-UVS-W/US (IP68 Submergence Cleaning)

* AK-C8100FA-UVS-W/US (upon request)

Keyboards for Medical, Industry
Commerce and Service

